

Italgin

✔ **ITALGIN** is a chromatic alginate for high precision impressions with high elasticity and resistance to tear.



✔ **Performance Highlights**

Dust free - Lead and cadmium free	Highest Compatibility
chromatic phase indicator	preparation of the impression is greatly simplified
Pleasant minty taste	For more patient comfort
Short intra-oral setting time	For real time-saving

✔ **Technical Characteristic**

Mixing time	Working time	Initial setting time	Time in mouth	Reproduction of details	Elastic recovery
60"	1'50"	1'55"	Minimum of 60"	50 µm	95.5%

Using deionized water at 23°C

Italgin



✓ **Bag Net Weight 450 gram**

✓ **Flavor: Mint**

✓ **Indications**

- Study models
- Fixed and removable prothesis
- Temporary impressions

✓ **Warning**

- Wear a mask, safety glasses, gloves and protective clothing. While being “Dust Free”, take care not to inhale the powder. Close the bag immediately after use. Don't use on patients allergic to MINT

✓ **Usage and storage tips**

- During the preparation of the impression, pour first water then powder inside the bowl. This will reduce the air bubbles and improve the mixture quality. To mix: find BMS spatula and Bowl at p. 76.
- During the filling phase of the tray make sure that the material comes out of the perforations, in this way the grip between the material and the tray will be stronger.
- Rinse with water and disinfect impression before storing.
- Store in as dry place and not exposed to the sunlight, keep Italgin in the original envelope close with a clip.



17,5 ml	+	9 g.
1 x		1 x
35 ml	+	18 g.
2 x		2 x

***Measurements taken at 23° C with deionized water.**

(*) Mixing time 60" (*) Working time 1'50" (*) Initial setting time 1'55" (+ 60" minimum time in mouth)