

kriovit

ICE SPRAY

- ✓ suitable for vitality testing of teeth and for quick cooling impression/casting materials i.e. hydrocolloid etc
- ⚠ **Do not use directly on gingival mucous, danger of freezing.**

✓ Main characteristics

- With pleasant scent of orange
- Precise spraying using application tubes

✓ Instructions for use

spray the object to be cooled from distance of approximately 2-3 cm for 1-2 seconds

✓ Composition

Propane
Butane pure
Ethanol
Orange juice oil

✓ Packing

200 mL

✓ Storage

- Keep away from heat, hot surfaces, sparks, open flames and other ignition sources.
- No smoking.
- P410+P412: Protect from sunlight.
- Do not expose to temperatures exceeding 50°C/122°F.
- Keep out of reach of children

