

EXPOSÉ

Disposable Caries Indicator

✓ EXPOSÉ IS A CARIES INDICATOR ON A STICK!

Our indicator solution is pre-dried right onto our Benda[®] Micro applicator. When applied to the tooth, Exposé stains the infected dentin. Simply touch Exposé on a wet tooth (or wet with a 3-way syringe) and apply it to the suspected lesion and rinse. If the lesion is stained, it is carious. Exposé will not stain healthy dentin or enamel.

✓ WITH EXPOSÉ YOU CAN BE SURE

- Ideal for verifying 'sticks' (cavitated or non-cavitated)
- Selective caries removal
- Ensure all carious material is removed from a prep or wall of a canal during a core build-up
- Exposé helps ensure that healthy dentin completely surrounds any remaining carious material

✓ FEATURES

- Pre-dried indicator solution so there is no inadvertent staining
- Two head sizes to make access easy – any tooth, any location, any patient
- Bendable head stays where you want it for even easier access
- Safe, patient-friendly FD&C dye in a propylene glycol base
- A cost so low that Exposé should be on every patient tray, every time

✓ Size

Fine (1.5 mm point)

Extra-Fine (1.0 mm point)

✓ Easier to see decay on dentin

✓ Differentiate between decay and pulp in deep preps

✓ Clearly see decay on preps where dentin decay is tooth-colored

