

CEMBEST

Glass Ionomer Luting Cement

✓ Characteristics and Advantages

- Good mixing properties: easy to mix (hand mix) resulting in a creamy consistency liner.
- Good pulp and tissues acceptance
- Excellent adhesion: adheres to dentine and enamel without the use of chemical bonding liners
- Low film thickness & Radiopaque
- Fluoride release: controlled release of fluoride ions - which in turn reduce the number of harmful microorganisms in the buccal cavity and reduce cavities
- Completely saliva resistant (no need the use of rubber dams saving time and effort)
- An increased working time that, while still retaining a rapid set, it allows the nurse/dentist to have more working time to use the cement and to prepare the crown/bridge which requires luting, if necessary
- A high level of translucency (translucency develops in the mouth)

✓ Indications

Crowns, bridges, inlays, onlays, posts, orthodontic bands

✓ Technical Characteristic

Working Time (including Mixing)	Setting Time (Including Working)
3:30 min	7:30 min

✓ Packing

35 gram Powder
20 ml Liquid
Powder Scoop, Mixing Pad

